

# PurpleSelect™ CENTRIFUGE TUBES

- Available in two volumes – 15ml and 50ml
- Tamper proof graduation marks
- Compatible with all the standard centrifuge rotors
- Frosted white area for writing
- Smooth and easy to open and close cap
- DNase, RNase & Endotoxin free
- Two packing options available -bulk packed and individually packed



- **Latex-free**  
No Latex was used in the manufacturing of the product. This also includes all the packaging and shipping materials used in the production of the items.

## Compliance

- **Manufacturing**  
These tubes are manufactured in ISO 9001:2015 certified facility.
- **Raw material**  
The resins used for making the tubes are USP Class-VI certified, RoHS, FDA regulation CFR 21, and are free of Substances of Very High Concern (SVHC), REACH compliant, and are FDA approved for food contact. They are free from Bisphenol A (BPA), phthalates and cytotoxic effects.
- **Sterilization**  
Product has been sterilized by **gamma** irradiation & released per the requirements of ISO 11137 - "Sterilization of healthcare products - radiation". Products meet a minimum Sterility Assurance Level (SAL) of  $10^{-6}$ .
- **TSE-BSE free**  
The product is free from and manufactured from materials that do not contain any raw materials or substances derived from animal origin as defined in EC Directive 1069/2009/EC and EC Directive 142/2011/EC.

## Technical Specifications

Specification	15ml tubes	50ml tubes	50ml self-standing tubes
Product codes	AVL001, AVL001I	AVL002, AVL002I	AVL002SS
Maximum volume	15ml	50ml	50ml
Graduation	starts from 1.5ml at every 0.5ml	starts from 5ml at every 2.5ml	starts from 5ml at every 2.5ml
MOC of tube	Virgin Polypropylene	Virgin Polypropylene	Virgin Polypropylene
MOC of cap	HDPE	HDPE	HDPE
Max. G-force	17000g	20000g	-
Maximum temperature	-80 to 121°C, 121°C for 15 min at 15psi		
Tube height (w/o cap)	119.30mm	114.50mm	114.50mm
Tube diameter (w/o cap)	17.20mm	29.80mm	29.80mm
Cap height	10mm	13mm	13mm
Cap diameter	21.80mm	35mm	35mm
Writing area	Alcohol resistant & tamper proof		

## Ordering Information

Name	Catalog No	Pack size
Centrifuge Tubes, 15ml, <b>Sterile</b> With screw cap and conical bottom DNAse-RNAse-Endotoxin-free Autoclavable, <b>Bulk-packed</b> , Gamma irradiated	AVL001-20X25NO	25 tubes per pack (500 tubes per case)
Centrifuge Tubes, 15ml, <b>Sterile</b> With screw cap and conical bottom DNAse-RNAse-Endotoxin-free Autoclavable, <b>Individually-packed</b> , Gamma irradiated	AVL001I-20X25NO	25 tubes per pack (500 tubes per case)
Centrifuge Tubes, 50ml, <b>Sterile</b> With screw cap and conical bottom DNAse-RNAse-Endotoxin-free Autoclavable, <b>Bulk-packed</b> , Gamma irradiated	AVL002-20X25NO	25 tubes per pack (500 tubes per case)
Centrifuge Tubes, 50ml, <b>Sterile</b> With screw cap and conical bottom DNAse-RNAse-Endotoxin-free Autoclavable, <b>Individually-packed</b> , Gamma irradiated	AVL002I-10X25NO	25 tubes per pack (250 tubes per case)
Centrifuge Tubes, 50ml, <b>Sterile</b> With screw cap and <b>Self-standing</b> DNAse-RNAse-Endotoxin-free Autoclavable, <b>Bulk-packed</b> , Gamma irradiated	AVL002SS-20X25NO	25 tubes per pack (500 tubes per case)



Discover More Products to Simplify Your Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at [info@purplebiolabs.com](mailto:info@purplebiolabs.com)/[support@purplebiolabs.com](mailto:support@purplebiolabs.com)



11/25\_00

Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.