

PurpleSelect™ CENTRIFUGE TUBES (AMBER COLOR TUBES)

- Amber color tubes for light sensitive samples
- Available in two volumes – 15ml and 50ml
- Tamper proof graduation marks
- Compatible with all the standard centrifuge rotors
- Frosted white area for writing
- Smooth and easy to open and close cap
- DNase, RNase & Endotoxin free
- Available in bulk packing

Compliance

- **Manufacturing**
These tubes are manufactured in ISO 9001:2015 certified facility.
- **Raw material**
The resins used for making the tubes are USP Class-VI certified, RoHS, FDA regulation CFR 21, and are free of Substances of Very High Concern (SVHC), REACH compliant, and are FDA approved for food contact. They are free from Bisphenol A (BPA), phthalates and cytotoxic effects.
- **Sterilization**
Product has been sterilized by **gamma** irradiation & released per the requirements of ISO 11137 - "Sterilization of healthcare products - radiation". Products meet a minimum Sterility Assurance Level (SAL) of 10^{-6} .
- **Latex-free**
No Latex was used in the manufacturing of the product. This also includes all the packaging and shipping materials used in the production of the items.

Technical Specifications

| Specification | 15ml tubes | 50ml tubes | 50ml self-standing tubes | 50ml Sterile, Self-standing tubes |
|---------------------|---|----------------------|--------------------------|-----------------------------------|
| Product codes | PBL023N | PBL024N | PBL025NS | PBL026SS |
| Maximum volume | 15ml | 50ml | 50ml | 50ml |
| MOC of tube | Virgin Polypropylene | Virgin Polypropylene | Virgin Polypropylene | Virgin Polypropylene |
| MOC of cap | HDPE | HDPE | HDPE | HDPE |
| Max. G-force | 17000g | 20000g | - | - |
| Maximum temperature | -80 to 121°C, 121°C for 15 min at 15psi Alcohol resistant & tamper proof | | | |
| Writing area | | | | |

Ordering Information

| Name | Catalog No | Pack size |
|---|------------|--|
| Centrifuge Tubes, 15ml, With screw cap and amber color DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated | PBL023N | 25 tubes per pack 20 packs per case Total 500 tubes per case |
| Centrifuge Tubes, 50ml, With screw cap and amber color DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated | PBL024N | 25 tubes per pack 20 packs per case Total 500 tubes per case |
| Centrifuge Tubes, 50ml self standing, With screw cap and amber color DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated | PBL025NS | 25 tubes per pack 20 packs per case Total 500 tubes per case |
| Centrifuge Tubes, 50ml self standing, Sterile With screw cap and amber color DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated | PBL026SS | 25 tubes per pack 20 packs per case Total 500 tubes per case |



Discover More Products to Simplify Your Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at info@purplebiolabs.com/support@purplebiolabs.com



11/25_00

Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.