

PurpleSelect™ CRYOVIALS

- Available in volumes – 1ml, 1.8ml, 3.5ml and 4.5ml
- V-shaped base ensures better sample collection
- Big labelling area for clear and easy sample marking
- Fast-lock cap for secure, leak-free sealing.
- Smooth and easy to open and close cap
- Stable hepta-foot design
- External threading design to prevent contamination
- DNase & RNase & Endotoxin free



Compliance

- **Manufacturing**
These tubes are manufactured in ISO 9001:2015 certified facility.
- **Raw material**
 - a) Meets strict US FDA safety and quality standards.
 - b) Safe for use in lab, medical, and life-science applications.
 - c) Ensures the material is non-toxic and free from harmful contaminants.



- **Latex-free**
No Latex was used in the manufacturing of the product. This also includes all the packaging and shipping materials used in the production of the items.

Technical Specifications

Specification	1ml cryovial	1.8ml cryovial	3.5ml cryovial	4.5ml cryovial
Product codes	PBL019SB	PBL020SB	PBL021SB	PBL022SB
Maximum volume	1ml	1.8ml	3.5ml	4.5ml
MOC of tube	Virgin Polypropylene	Virgin Polypropylene	Virgin Polypropylene	Virgin Polypropylene
Sterile/Non-Sterile	Sterile	Sterile	Sterile	Sterile
Maximum temperature	-196°C to 121°C, 121°C for 15 min at 15psi			
Writing area	Alcohol resistant & tamper proof			

Ordering Information

Name	Catalog No	Pack size
Cryovials, 1ml, Sterile With screw cap and hepta foot DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated	PBL019SB	100 vials per pack, 10 packs per case (Total 1000 vials per case)
Cryovials, 1.8ml, Sterile With screw cap and hepta foot DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated	PBL020SB	100 vials per pack, 10 packs per case (Total 1000 vials per case)
Cryovials, 3.5ml, Sterile With screw cap and hepta foot DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated	PBL021SB	100 vials per pack, 10 packs per case (Total 1000 vials per case)
Cryovials, 4.5ml, Sterile With screw cap and hepta foot DNAse-RNAse-Endotoxin-free Autoclavable, Bulk-packed , Gamma irradiated	PBL022SB	100 vials per pack, 10 packs per case (Total 1000 vials per case)



Discover More Products to Simplify Your Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at info@purplebiolabs.com/support@purplebiolabs.com



11/25_00

Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.