

PurpleSelect™ PASTEUR PIPETTES

- Available in volumes – 1ml and 3ml
- Built from durable, medical-grade LDPE for safe handling
- Soft, easy-squeeze bulb ensures smooth pipetting control
- Polished opening provides clean, accurate liquid release
- Expertly produced for reliable and efficient performance
- Clearly marked graduations support precise volume reading
- Fine molding quality gives a sleek, ultra-smooth inner surface
- DNase, RNase & Endotoxin free



Compliance

Manufacturing

- These pipettes are manufactured under ISO Class 8 manufacturing facility according to ISO 14644-1
- The material used is high-purity, virgin-grade polypropylene that contains absolutely no metal traces, ensuring outstanding clarity, reliable performance, and complete safety for sensitive laboratory applications.

Technical Specifications

Specification	Pasteur Pipette	Pasteur Pipette	Pasteur Pipette	Pasteur Pipette
Product codes	PBL017B	PBL017S	PBL018B	PBL018S
Maximum volume	1ml	1ml	3ml	3ml
Sterile/Non-Sterile	Non-Sterile	Sterile	Non-Sterile	Sterile
Packing	Bulk	Individually wrapped	Bulk	Individually wrapped
MOC of tube	Grade Virgin Polypropylene	Grade Virgin Polypropylene	Grade Virgin Polypropylene	Grade Virgin Polypropylene

Ordering Information

Name	Catalog No	Pack size
Pasteur Pipettes, 1ml, Bulk Transparent & clear DNAse-RNAse-Endotoxin-free	PBL017B	500 pcs per pack, 10 packs per case (Total 5000 pcs per case)
Pasteur Pipettes, 1ml, Individually Wrap Sterile Transparent & clear DNAse-RNAse-Endotoxin-free	PBL017S	300x1 per pack, 10 packs per case (Total 3000 pcs per case)
Pasteur Pipettes, 3ml, Bulk Transparent & clear DNAse-RNAse-Endotoxin-free	PBL018B	500 pcs per pack, 10 packs per case (Total 5000 pcs per case)
Pasteur Pipettes, 3ml, Individually Wrap Sterile Transparent & clear DNAse-RNAse-Endotoxin-free	PBL018S	300x1 per pack, 10 packs per case (Total 3000 pcs per case)



Discover More Products to Simplify Your Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at info@purplebiolabs.com/support@purplebiolabs.com

11/25_00



Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.