

PurpleSelect™ SEROLOGY PIPETTES

- Available in volumes – 1ml, 2ml, 5ml and 10ml
- **Accuracy**
 - Delivers precise volume control within $\pm 2\%$.
- **PET Filter Plug**
 - PET filter acts as a barrier to protect samples and avoid contamination.
- **Individual Packaging**
 - Individually sealed in traceable peel-back packs for easy verification with visible lot details.
- **Bidirectional Graduations**
 - Graduations readable from both sides for flexible handling.
- **Color Coding**
 - Color bands simplify size recognition.
- DNase, RNase, Endotoxin free



Compliance

- **Manufacturing**

These tubes are manufactured in ISO 9001:2015 and ISO 13485:2016 certified facility.
- **Raw material**

The resins used for making the tubes are USP Class-VI certified, RoHS, FDA regulation CFR 21, and are free of Substances of Very High Concern (SVHC), REACH compliant, and are FDA approved for food contact. They are free from Bisphenol A (BPA), phthalates and cytotoxic effects.
- **Sterilization**

Product has been sterilized by **gamma** irradiation & released per the requirements of ISO 11137 - "Sterilization of healthcare products - radiation". Products meet a minimum Sterility Assurance Level (SAL) of 10^{-6} .
- **Latex-free**

No Latex was used in the manufacturing of the product. This also includes all the packaging and shipping materials used in the production of the items.

Technical Specification

Specification	1ml cryovial	2.0ml cryovial	5.0ml cryovial	10ml cryovial
Product codes	PBL013S	PBL014S	PBL015S	PBL016S
Maximum volume	1ml	2ml	5ml	10ml
MOC of tube	Virgin Polystyrene	Virgin Polystyrene	Virgin Polystyrene	Virgin Polystyrene
Packing	Individually wrapped	Individually wrapped	Individually wrapped	Individually wrapped
Color Code	Yellow	Green	Blue	Orange

Ordering Information

Name	Catalog No	Pack size
Serological Pipettes, 1ml, Sterile DNAse-RNAse-Endotoxin-free Individually wrapped, Gamma irradiated	PBL013S	100 pcs per pack, 6 packs per case (Total 600 pcs per case)
Serological Pipettes, 2ml, Sterile DNAse-RNAse-Endotoxin-free Individually wrapped, Gamma irradiated	PBL014S	100 pcs per pack, 6 packs per case (Total 600 pcs per case)
Serological Pipettes, 5ml, Sterile DNAse-RNAse-Endotoxin-free Individually wrapped, Gamma irradiated	PBL015S	100 pcs per pack, 6 packs per case (Total 600 pcs per case)
Serological Pipettes, 10ml, Sterile DNAse-RNAse-Endotoxin-free Individually wrapped, Gamma irradiated	PBL016S	100 pcs per pack, 6 packs per case (Total 600 pcs per case)



Discover More Products to Simplify Your
Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at
info@purplebiolabs.com/support@purplebiolabs.com

11/25_00



Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.