

PurpleSelect™ CELL STRAINERS

- Made up of medical grade polypropylene frame and imported nylon mesh
- Available in three pore sizes – 40, 70 and 100µm.
- Colour coding for each pore size for easy identification
- Handle on one side for easy operation
- Compatible with most of the brands of 50ml tubes
- Gamma sterilized and individually packed
- Endotoxin free



Ordering Information

Name	Catalog No	Pack size
Cell strainers, Pore size: 70 µm Sterile, Individually packed	PBL008-50X1NO	50 individually packed cell strainers per pack
Cell strainers, Pore size: 40 µm Sterile, Individually packed	PBL009-50X1NO	50 individually packed cell strainers per pack
Cell strainers, Pore size: 100 µm Sterile, Individually packed	PBL010-50X1NO	50 individually packed cell strainers per pack



Discover More Products to Simplify Your
Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at
info@purplebiolabs.com/support@purplebiolabs.com



11/25_00

Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.