

PurpleCore™

Manual Petri Dish Revolving Table

- Hand-operated mechanism
- Heavy-duty rotating base for extended spinning
- Usable on both sides
- Suitable for Petri dishes up to 150 mm in diameter
- Rubberized anti-slip surface compatible with glass and plastic dishes



Description:

Designed for evenly spreading agar on Petri dishes, this device is commonly used in molecular biology and microbiology laboratories. The non-slip working surface allows for smooth and simple operation.

Operating Instructions:

Place the rotating platform on a level bench. Dispense the sample into the Petri dish, then rotate the plate while spreading the material evenly using a spreader.

Technical Specifications:

Specification	
Product Code	PBI001
Petri dish diameter (Φ)	60~150mm
Operation mode	Manual
Workplate	2 pcs
Weight	2.4 kg



Discover More Products to Simplify Your Laboratory Needs

Visit us at <https://purplebiolabs.com/>

For further queries, contact us at

info@purplebiolabs.com/support@purplebiolabs.com



Prior to usage, users are responsible for ensuring the suitability of our products for their intended application. Our products adhere solely to the information provided in this and other related Purple Biolabs publications. The content in this publication is derived from our research and development efforts and is, to the best of our knowledge, accurate and reliable. Purple Biolabs retains the right to modify specifications and product-related information at any time. Our products are designed for laboratory and research purposes only, unless otherwise specified. The statements herein should not be construed as a warranty, whether expressed or implied, and we do not accept liability for any patent violations.